

TSH FIA Rapid Test (Serum/Plasma) FIA-TSH-001

*A rapid test for detecting TSH in serum or plasma with the use of the
Biopanda Fluorescence Immunoassay Analyser.
For professional in vitro diagnostic use only.*

INTENDED USE

The Biopanda TSH FIA Rapid Test is based on Fluorescence Immunoassay for the quantitative determination of Thyroid Stimulating Hormone (TSH) in serum or plasma. Measurement of TSH is useful to aid in the screening of adult populations for primary hypothyroidism by medical professionals. It could also be used in screening neonates for hypothyroidism.

SUMMARY

Thyroid-stimulating hormone (also known as thyrotropin, thyrotropic hormone, TSH, or hTSH for human TSH) is a pituitary hormone that stimulates the thyroid gland to produce thyroxine(T₄), and then triiodothyronine(T₃) which stimulates the metabolism of almost every tissue in the body.[1] It is a glycoprotein hormone synthesized and secreted by thyrotrope cells in the anterior pituitary gland, which regulates the endocrine function of the thyroid.[2][3] TSH (with a half life of about an hour) stimulates the thyroid gland to secrete the hormone thyroxine(T₄), which has only a slight effect on metabolism. T₄ is converted to triiodothyronine(T₃), which is the active hormone that stimulates metabolism. About 80% of this conversion is in the liver and other organs, and 20% in the thyroid itself.[1] Laboratory testing of thyroid stimulating hormone levels in the blood is considered the best initial test for hypothyroidism.[4] It is important to note the statement from the Subclinical Thyroid Disease Consensus Panel: "There is no single level of serum TSH at which clinical action is always either indicated or contraindicated. The higher the TSH, the more compelling is the rationale for treatment. It is important to consider the individual clinical context (e.g. pregnancy, lipid profile, ATPO antibodies),"[5]

PRINCIPLE

The TSH FIA Rapid Test Cassette detects TSH based on Fluorescence Immunoassay. The sample moves through the strip from sample pad to absorbent pad. If the specimen contains TSH, it attaches to the TSH antibody which is conjugated with fluorescent microspheres. Then the complex will be captured by the capture antibody coated on the nitrocellulose membrane. The concentration of TSH in the sample correlates linearly with the fluorescence signal intensity captured on the T line. According to the fluorescence intensity of the test and product standard curve, the concentration of TSH in the sample can be calculated by the Biopanda Fluorescence Immunoassay Analyser to show TSH concentration in specimen.

REAGENTS

The test kit includes TSH antibody coated particles and TSH antibody coated on the membrane.

PRECAUTIONS

1. For professional *in vitro* diagnostic use only.
2. Do not use after the expiration date indicated on the package. Do not use the test if the foil pouch is damaged. Do not reuse.
3. This test contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not completely guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious, and handled observing usual safety precautions (e.g., do not ingest or inhale).
4. Avoid cross-contamination of specimens by using a new specimen collection container for each specimen obtained.
5. Do not eat, drink or smoke in the area where the specimens and tests are handled. Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout the procedure and follow standard procedures for proper disposal of specimens. Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are assayed.

6. Do not interchange or mix reagents from different lots.
7. Humidity and temperature can adversely affect results.
8. Used testing materials should be discarded in accordance with local regulations.
9. Read the entire procedure carefully prior to any testing.
10. The Biopanda TSH FIA Rapid Test should only be used with the Biopanda Fluorescence Immunoassay Analyser by approved medical professionals.

STORAGE AND STABILITY

1. The kit should be stored at 4-30°C until the expiry date printed on the sealed pouch.
2. The test must remain in the sealed pouch until use.
3. Do not freeze.
4. Care should be taken to protect the components of the kit from contamination. Do not use if there is evidence of microbial contamination or precipitation. Biological contamination of dispensing equipment, containers or reagents can lead to false results.

KIT COMPONENTS

- 10 x foil wrapped TSH test cassettes
- 10 x Specimen collection tubes with dilution buffer
- 1 x ID card (TSH)
- Package Insert

MATERIALS REQUIRED BUT NOT PROVIDED

- Timer
- Centrifuge
- Pipette
- Biopanda Fluorescence Immunoassay Analyser

SPECIMEN COLLECTION AND PREPARATION

PREPARATION

1. Before performing the test, please make sure that all components are brought to room temperature (15-30°C). Cold buffer solution or moisture condensation on the membrane can lead to invalid test results.
2. Take a tube with buffer solution out of the kit. Document patients name or ID on it. Open the screw cap.

BLOOD SAMPLE TAKING

1. Collect the specimen according to standard procedures.
2. Do not leave specimens at room temperature for prolonged periods. Serum and plasma specimens may be stored at 2-8 °C for up to 1 day, for long term storage, specimens should be kept below -20 °C.
3. Bring specimens to room temperature prior to testing. Frozen specimens must be completely thawed and mixed well prior to testing. Avoid repeated freezing and thawing of specimens.
4. EDTA and Heparin sodium, can be used as the anticoagulant for collecting the specimen.

SAMPLE DILUTION / SAMPLE STABILITY

1. Transfer **75 µL of serum or plasma** to the buffer tube with a micro pipette.
2. Close the tube and shake the sample by hand vigorously for approximately **10 seconds** to mix the sample and dilution buffer.
3. Let the diluted sample homogenize for approximately 1 minute.
4. The diluted sample can then be used immediately or stored for up to 8 hours.

DIRECTIONS FOR USE

Refer to the Biopanda Fluorescence Immunoassay Analyser Operation Manual for the complete instructions on use of the Test. The test should be conducted at room temperature.

Allow the test, specimen, buffer and/or controls to reach room temperature (15-30°C) prior to testing.

1. Turn on the Analyser. Then according to the user requirement, select "Standard test" or "Quick test" mode.
2. Take out the ID card and insert it into the Analyser port.
3. **Serum/plasma:** Transfer 75 µl of serum/plasma into the buffer tube, mix the specimen and the buffer well.
4. **Add diluted specimen with a Pipette:** Pipette 75 µl of diluted specimen into the sample well of the test cassette. Start the timer at the same

time.

- There are two test modes for the Biopanda Fluorescence Immunoassay Analyser; Standard Test mode and Quick Test mode. Please refer to the user manual of the Biopanda Fluorescence Immunoassay Analyser for details.

“Quick test” mode: Insert the test cassette into the Analyser at 15 minutes after sample application and click "New Test", the Analyser will automatically give the test result after a few seconds.

“Standard test” mode: Insert the test cassette into the Analyser immediately after sample application, click "New test" at the same time, the Analyser will automatically count down the 15 minutes. After the countdown, the Analyser will give the result at once.

INTERPRETATION OF RESULTS

The result of tests for TSH is calculated by the Biopanda Fluorescence Immunoassay Analyser and displays the result on the screen. For additional information, please refer to the user manual of the Biopanda Fluorescence Immunoassay Analyser.

Linearity range of the Biopanda TSH FIA Rapid Test is 0.1-100 µIU/ml.

QUALITY CONTROL

Each Biopanda TSH FIA Rapid Test Cassette contains an internal control that satisfies routine quality control requirements. This internal control is performed each time a patient sample is tested. This control indicates that the test device was inserted and read properly by the Biopanda Fluorescence Immunoassay Analyser. An invalid result from the internal control causes an error message on the Biopanda Fluorescence Immunoassay Analyser indicating that the test should be repeated. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

LIMITATIONS

- The Biopanda TSH FIA Rapid Test Cassette is for professional *in vitro* diagnostic use, and should only be used for the quantitative detection of TSH.
- The Biopanda TSH FIA Rapid Test Cassette will only indicate the presence of TSH antigen in the specimen and should not be used as the sole criterion for evaluating thyroid function.
- As with all diagnostic tests, a confirmed diagnosis should only be made by a physician after all clinical and laboratory findings have been evaluated.
- The results of the Biopanda TSH FIA Rapid Tests are based on measuring the levels of TSH in a specimen. It should not be used as the sole criterion for treatment decisions. If the result is positive, other clinical findings and alternative test methods are recommended to reach proper medical treatments.

EXPECTED RESULTS

| Concentrations | Clinical Reference |
|----------------|--------------------|
| <20 µIU/ml | Normal Neonatal |
| <10 µIU/ml | Normal Children |
| <5 µIU/ml | Normal Adult |

PERFORMANCE CHARACTERISTICS

- ACCURACY:** The test deviation is $\leq \pm 15\%$.
- DETECTION RANGE:** 0.1-100 µIU/ml.
- LINEARITY RANGE:** 0.1-100 µIU/ml, $R \geq 0.990$
- PRECISION**

INTRA-LOT PRECISION

Within-run precision has been determined by using 10 replicates of 2 specimens containing 5.0 µIU/ml and 20 µIU/ml of TSH. C.V. is $\leq 15\%$.

INTER-LOT PRECISION

Between-run precision has been determined by using 10 replicates for each of three lots using 2 specimens containing 5.0 µIU/ml and 20 µIU/ml of TSH. C.V. is $\leq 15\%$.

REFERENCES

- Merck Manual of Diagnosis and Therapy, Thyroid gland disorders.
- The American Heritage Dictionary of the English Language, Fourth Edition. Houghton Mifflin Company. 2006. ISBN0-395-82517-2.

- Sacher R, Richard A. McPherson (2000). Widmann's Clinical Interpretation of Laboratory Tests, 11th ed. F.A. Davis Company. ISBN0-8036-0270-7.
- So, M; MacIsaac, RJ; Grossmann M (August 2012). "Hypothyroidism". Australian Family Physician 41 (8): 556–62.
- Surkset. al., JAMA 291:228, 2004.

Index of Symbols

| | | | | | |
|--|---|--|-----------------|--|------------------------------|
| | Manufacturer | | Tests per kit | | Do not reuse test |
| | <i>In vitro</i> diagnostic medical device | | Expiration date | | Catalogue number |
| | Store between 4-30°C | | Lot Number | | Consult instructions for use |
| | Do not use if package is damaged | | | | |

Thank you for purchasing Biopanda's TSH FIA Rapid Test. Please read this manual carefully before operating to ensure proper use.



Biopanda Reagents Ltd.

Unit 14 Carrowreagh Business Park
Carrowreagh Road
Belfast, BT16 1QQ
United Kingdom
Tel: +44 (0) 28 95438774
E-mail: info@biopanda.co.uk
Website: www.biopanda.co.uk

Effective date: 09/09/2021