

AMH FIA Rapid Test Device



Easy to use rapid test for the quantitative detection of Antimullerian hormone (AMH) in whole blood, serum, or plasma samples.

Measurement of AMH is used to aid in the assessment of ovarian reserve function, menopause, and polycystic ovary syndrome (PCOS).

Test Features:

- Rapid results in only 15 minutes
- Whole blood, serum, or plasma samples
- Detection range: 0.05-25 ng/ml
- Storage Temperature: 4-30°C
- Use with the Fluorescence Immunoassay Device (BR-FIA-2000)

Ordering information:

PRODUCT NAME	CATALOGUE No	KIT SIZE
AMH FIA Rapid Test Device	FI2-AMH-001	25 Tests

BioPanda Reagents Ltd.

Unit 14, Carrowreagh Business Park, Carrowreagh Road, Dundonald
Belfast, BT16 1QQ, Northern Ireland, United Kingdom
Tel : +44 (0)28 95438774 Email: info@biopanda.co.uk Web: www.biopanda.co.uk

