

β -HCG FIA Rapid Test Device



Easy to use rapid test for the quantitative detection of β -human chorionic gonadotropin (β -hCG) in whole blood, serum, or plasma samples.

Measurement of β -hCG is primarily used to aid in pregnancy detection but can also be used as a tumour marker in not pregnant individuals.

Test Features:

- Rapid results in only 15 minutes
- Whole blood, serum, or plasma samples
- Detection range: 2-200,000 mIU/ml
- Storage Temperature: 4-30°C
- Use with the Fluorescence Immunoassay Device (BR-FIA-2000)

Ordering information:

PRODUCT NAME	CATALOGUE No	KIT SIZE
β -HCG FIA Rapid Test Device	FI2-BHCG-001	25 Tests

BioPanda Reagents Ltd.

Unit 14, Carrowreagh Business Park, Carrowreagh Road, Dundonald
Belfast, BT16 1QQ, Northern Ireland, United Kingdom
Tel : +44 (0)28 95438774 Email: info@biopanda.co.uk Web: www.biopanda.co.uk

