

# CA19-9 FIA Rapid Test Device



Easy to use rapid test for the quantitative detection of CA19-9 in whole blood, serum, or plasma samples.

Measurement of CA19-9 is used as an aid in the diagnosis and treatment monitoring of pancreatic cancer and other gastrointestinal malignant tumours.

Please note: A confirmed diagnosis should only be made by a physician after all clinical and laboratory findings have been evaluated.

## Test Features:

- Rapid results in only 15 minutes
- Whole blood, serum, or plasma samples
- Detection range: 2-2000 U/mL
- Storage Temperature: 4-30°C
- Use with the Fluorescence Immunoassay Device (BR-FIA-2000)

## Ordering information:

PRODUCT NAME	CATALOGUE No	KIT SIZE
CA19-9 FIA Rapid Test Device	FI2-CA199-001	25 Tests

### BioPanda Reagents Ltd.

Unit 14, Carrowreagh Business Park, Carrowreagh Road, Dundonald  
Belfast, BT16 1QQ, Northern Ireland, United Kingdom  
Tel : +44 (0)28 95438774 Email: info@biopanda.co.uk Web: www.biopanda.co.uk

