

Calprotectin FIA Rapid Test Device



Easy to use rapid test for the quantitative detection of Calprotectin in human faecal samples.

Measurement of Calprotectin is used as an aid in the assessment of intestinal inflammation.

Test Features:

- Rapid results in only 15 minutes
- Human faecal samples
- Detection range: 1-1000 $\mu\text{g/g}$
- Storage Temperature: 4-30°C
- Use with the Fluorescence Immunoassay Device (BR-FIA-2000)

Ordering information:

PRODUCT NAME	CATALOGUE No	KIT SIZE
Calprotectin FIA Rapid Test Device	FI2-CAL-001	25 Tests

BioPanda Reagents Ltd.

Unit 14, Carrowreagh Business Park, Carrowreagh Road, Dundonald
Belfast, BT16 1QQ, Northern Ireland, United Kingdom
Tel : +44 (0)28 95438774 Email: info@biopanda.co.uk Web: www.biopanda.co.uk

