

FOB FIA Rapid Test Device



Easy to use rapid test for the quantitative detection of Faecal Occult Blood (FOB) in human faecal samples.

Measurement of FOB is used as an aid in the diagnosis of gastrointestinal problems.

Test Features:

- Rapid results in only 15 minutes
- Human faecal samples
- Quantitative results
- Storage Temperature: 4-30°C
- Use with the Fluorescence Immunoassay Device (BR-FIA-2000)

Ordering information:

| PRODUCT NAME | CATALOGUE No | KIT SIZE |
|---------------------------|--------------|----------|
| FOB FIA Rapid Test Device | FI2-FOB-001 | 25 Tests |

BioPanda Reagents Ltd.

Unit 14, Carrowreagh Business Park, Carrowreagh Road, Dundonald
Belfast, BT16 1QQ, Northern Ireland, United Kingdom
Tel : +44 (0)28 95438774 Email: info@biopanda.co.uk Web: www.biopanda.co.uk

