

H. Pylori Ag FIA Rapid Test Device



Easy to use rapid test for the quantitative detection of *H. pylori* antigen in human faecal samples.

Measurement of *H. pylori* antigen is used as an aid in the diagnosis of *H. pylori* infections.

Test Features:

- Rapid results in only 15 minutes
- Human faecal samples
- Detection range: 2-1500 ng/mL
- Storage Temperature: 4-30°C
- Use with the Fluorescence Immunoassay Device (BR-FIA-2000)

Ordering information:

PRODUCT NAME	CATALOGUE No	KIT SIZE
H. Pylori Ag FIA Rapid Test Device	FI2-HPAG-001	25 Tests

BioPanda Reagents Ltd.

Unit 14, Carrowreagh Business Park, Carrowreagh Road, Dundonald
Belfast, BT16 1QQ, Northern Ireland, United Kingdom
Tel : +44 (0)28 95438774 Email: info@biopanda.co.uk Web: www.biopanda.co.uk

