

Rotavirus FIA Rapid Test Device



Easy to use rapid test for the qualitative detection of Rotavirus antigen in human faecal samples.

Measurement of Rotavirus antigen is used as an aid in the diagnosis of Rotavirus infections

Test Features:

- Rapid results in only 15 minutes
- Human faecal samples
- Positive or negative results
- Storage Temperature: 4-30°C
- Use with the Fluorescence Immunoassay Device (BR-FIA-2000)

Ordering information:

| PRODUCT NAME | CATALOGUE No | KIT SIZE |
|---------------------------------|--------------|----------|
| Rotavirus FIA Rapid Test Device | FI2-RTV-001 | 25 Tests |

BioPanda Reagents Ltd.

Unit 14, Carrowreagh Business Park, Carrowreagh Road, Dundonald
Belfast, BT16 1QQ, Northern Ireland, United Kingdom
Tel : +44 (0)28 95438774 Email: info@biopanda.co.uk Web: www.biopanda.co.uk

